Synthesis of Aspirin

Synthesis Purification Characterization


Aspirin: some background

Patented by Bayer in 1893
One of the oldest drugs
One of the most consumed drugs (Production in the US is 10 million Kg/year)

Aspirin: Biological activity

- Analgesic (painkiller)
- Antipyretic (fever reducer)
- Anti-inflammatory (inhibition of the synthesis of prostaglandins)
- Side effects: gastric irritation, bleeding Apparition of new analgesics (Tylenol)

Aspirin: The Molecule



acetyl salicylic acid (aspirin)

Organic background Alcohols

Alcohols

- R—OH
- <u>Aliphatic alcohols</u>: Hydroxyl derivatives of saturated hydrocarbons


H₃C—OH

methanol




ethanol

Enols: Hydroxyl derivatives of unsaturated hydrocarbons




enol






Organic background Esterification



This is an equilibrium!

Organic background Esterification mechanism




Organic background Acidity: a review

<u>Bronsted-Lowry acid</u>
 Substance that can donate a proton
 <u>Lewis acid</u>
 Substance that accepts lone-pair electrons


• $AH + H_2O \rightarrow A^- + H_3O^+$ $K_a = [H_3O^+][A^-] / [AH]$ $pK_a = -Log K_a$

Acid	pK _a
water	14.0
Acetic acid	4.75
Benzoic acid	4.19
Formic acid	3.75
Sulfuric acid	1

Reaction



Reaction mechanism



Safety

 Consult MSDS for safety information on all the chemicals you plan to use.

■ MSDS are available from the 2507 webpage.

Synthesis procedure

- Mix salicylic acid (solid) with acetic anhydride (liquid) and a catalytic amount of H₂SO₄.
- Heat for 10 min at 90° C.
- Cool down to room temperature and add water (crystallization solvent).
- Induce crystallization.
- **Filter the crystals.**

Purification by Recrystallization WHAT IS

RECRYSTALLIZATION?

- Rapid and convenient way of purifying a solid organic compound
- > The material to be purified is dissolved in the hot appropriate solvent
- > As the solvent cools, the solution become saturated with respect to the substance, which then crystallize
- Impurities stay in solution

METHOD

- Choose the solvent
- > Dissolve the solute
- (Filter suspended solids)
- Crystallize the solute
- Collect and wash the crystals
- Dry the crystals

Assessment of purity: Looking for properties

Physical properties

Chemical properties

- Physical appearance
- Solubility
- Melting point
- Infrared spectra

- Percent yield
- Potentiometric titration
- Phenol test

Assessment of purity: Physical appearance

- Physical state (solid? liquid? ..)Color
- Odor
- Texture
- Homogeneity in composition

Assessment of purity: Solubility

- In a test tube, transfer a small amount of your product (end of a spatula)
- Add ~1/2mL of solvent
- Determine solubility at room temperature
- If not, gently heat and determine solubility at higher temperature
- Test several organic solvents and compare with litterature

Assessment of purity: Melting point

- Melting point characterize the compound
 Range of the melting point indicate the purity of the compound
- Method: scanning of the temperature until melting occurs. Determine start and end temperature of melting.
- Don't heat too fast!! (<1^oC / min)

Assessment of purity: Infrared Spectra

- Dissolve the product in chloroform and obtain the IR spectrum.
- Assign the peaks

- Do you see impurities?
- Compare with the reference spectrum.

Assessment of purity: Percent yield

- Definition: % Yield = n_{exp}/n_{th}
- Write down the reaction with the correct stoechiometry. Ex: $A+B \rightarrow C$
- Prepare a chart like below and calculate yield.

m _A (g) or V _A (mL)	n _A (mol)	m _B (g) or V _B (mL)	n _B (mol)	m _C (g) theoretical	n _C (mol) theoretical	m _c (g) experimental	n _C (mol) experimental

Assessment of purity: potentiometric titration

- <u>Objective</u>: Determine the pKa at the half equivalence point.
- Procedure: Titrate potentiometrically the acetylsalicylic acid with a strong base and find pKa



Assessment of purity: Phenol test and visible spectroscopy

- Phenol react with $FeCl_3$ (aq) to give a deep purple complex.
- Phenol is not present in the product but in one of the reactant.
- This test indicate the presence of **unreacted starting material** (quantitative analysis possible via visible spectroscopy).

