

8—Fixed Space Maintainer


Unilateral band-and-loop space maintainers are passive orthodontic appliances designed to preserve the space created by premature loss or extraction of teeth. Variations of this appliance include the use of a preformed stainless steel crown instead of a band or the use of a subgingival shoe in order to guide an unerupted tooth into proper alignment.

Received 8/24/90

DENTAL SCHOOL UNIVERSITY OF MARYLAND

Student DR. SMITH Date 8/23/90

Patient CARA WARFIELD PR No. 9/8/90



Please construct a fixed maxillary band-and-loop space maintainer in .032-inch wire at tooth #3 area. Add occlusal rest (.025) as indicated to prevent vertical movement.

Thank you!

DISINFECTED

Instructor N. H. Smith Cubicle No. 4



SUBURBAN
ORTHODONTIC
LABORATORY, INC.
49 WEST NYACK ROAD
NANUET, NEW YORK 10954
914-623-1240 • 914-623-1241
1-800-245-7668
FAX: 914-623-3329



CERTIFIED DENTAL
LABORATORY

Doctor: _____

Address: _____ Phone: () _____

City: _____ State: _____ Zip: _____


Patient's Name: P, ALVAREZ



Color

SPECIAL INSTRUCTIONS:

Please fabricate
a band + loop
with band on
A



Signature: ABC License # 12345

Date Shipped: _____

Date Needed: _____

White: Original Yellow: Work Order Pink: Dr. Copy

PROVIDE TO THE LAB FOR ALL CASES

1. Bands on the model if bands required.
2. Model labeled with patient's name and doctor's name.
3. Appliance design drawn on the model.
4. Any teeth to be extracted removed from model by doctor to simulate healed tissues.
5. Send the white (original) and yellow (work order) forms to the lab along with the model. The Pink for is the Doctor's copy and should be inserted in the chart.